

CareSens® Dual Fact Sheet



An innovative, ergonomically designed meter which tests both blood glucose and ketone levels and exceeds ISO 15197:2013 standards of accuracy¹



- Reads both blood glucose and ketone test strips
- Bluetooth connectivity
- Exceeds ISO 15197:2013¹
- Auto coding minimises errors
- Stores 1,000 readings
- Only 0.4 µl sample size
- 5 second test time
- Backlit LCD screen
- Strip ejection function
- Easy to hold and grip

Blood Glucose Test Strips **£9.95**

1 Meter 2 strips

Ketone Test Strips **£9.95**

| | | | |
|-----------------------------------|-------------------------------|------------------------------------|---|
| Size 108 x 58 x 17 (mm) | Temperature 5-45 °C | Relative humidity 10-90% | Weight 72.6g (with batteries) |
|-----------------------------------|-------------------------------|------------------------------------|---|

| Item | Acceptance Criteria | Result | | |
|---|---|------------------------------|---------------------|---------------------|
| Glucose Concentration <5.55mmol/L | 95% of the measured glucose values falls within ±0.83mmol/L | within ±0.83 mmol/L | within ±0.55 mmol/L | within ±0.28 mmol/L |
| | | 99.5% | 96.6% | 66.7% |
| Glucose Concentration ≥5.55mmol/L | 95% of the measured glucose values falls within ±15% | within ±15% | within ±10% | within ±1% |
| | | 99.7% | 95.7% | 66.7% |
| 99% of glucose values fall within zones A and B of the Concensus Error Grid | | Within zones A and B 100% | | |

References:

1. Data on file, i-sens ISO 15197:2013 testing results, 2016



24/7 NON-AUTOMATED FREE SUPPORT

0800 881 5423

Spirit House, Saffron Way, Off Saffron Lane, Leicester, LE2 6UP

T: 0116 286 5000 E: info@spirit-healthcare.co.uk www.spirit-healthcare.co.uk